



Body part specific settings



Track your results and chart your progress



Deploy multiple iTENS devices at a time



Muscle massage and muscle rehabilitation modes

The itens™ is the world's first true wireless TENS device that works with your iPhone or Android device.



iTENS sends gentle electrical impulses through the wireless electrodes attached to pain sites



These tiny electrical signals elicit a sensory stimulation to the body's nervous system



Overriding pain signals from reaching the brain



Releases the body's natural pain fighting chemicals called endorphins



Tracks and charts the users pain reduction delta



iTENS

wearable electrotherapy



Visit itens.com for additional information, instructions, and additional details. info@itens.com

twitter.com/itensdevice

pinterest.com/itensdevice

facebook.com/itensdevice

Made for iPod iPhone iPad

Works with Android™

FREE download on app store

pain relief with the push of a button

The iTens™ is an electrotherapy device that merges modern technology with the proven results of “TENS therapy” to provide effective and lasting pain relief. The iTens harnesses the power of your smart phone, with an iPhone or Android app that wirelessly syncs with the Bluetooth based electrotherapy device.

- Fast acting + long lasting pain relief
- Syncs with smart phone app
- Lithium-ion battery provides 24-hours of use on a single charge
- FDA cleared with medical grade strength
- Thin, flexible + discreet for easy use under clothing
- No wires or bulky devices
- Reusable and replaceable gel pads for multiple applications
- No prescription necessary drug free with no side effects

Three easy steps

- 

1 PEEL N' STICK GEL
- 

2 PLACE ON PAIN
- 

3 LAUNCH APP

DEVICE KITS




SMALL WINGS
for ankle, knee, elbow, wrist and other small treatment areas

AVAILABLE IN:
ORANGE
#itensORGSM



LARGE WINGS
for back, shoulder and large treatment areas

AVAILABLE IN:
BLUE
#itensbLL



ACCESSORIES



SIZE	QUANTITY	PART #
SMALL WINGS 2" x 6.25"	3 pairs/Box	ITRODESM
LARGE WINGS 2.5" x 6.75"	3 pairs/Box	ITRODELG
LONG STRIP 1.5" x 7"	3 pairs/Box	ITRODESTR



LARGE WING

AVAILABLE IN:
BLUE
#itensBLLA



SMALL WING

AVAILABLE IN:
ORANGE
#itensORGsMA



LONG STRIP
for the sciatic nerve, tendons and other long treatment areas

AVAILABLE IN:
GREEN
#itensgrls



- Charger
- Charging cord